Isotopengenerator Cs-137/Ba-137m

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Isotope Generator Cs-137/Ba-137m

Handling Instructions

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■ Description and Application

The Eckert & Ziegler Nuclitec Cs-137/Ba-137m Isotope Generator has been specially designed for use in experiments in schools and universities to demonstrate the physical properties of radioactive decay.

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Isotope Generator Cs-137/Ba-137m

Handling Instructions

The Isotope Generator can produce up to 1,000 small aliquots of a short-lived radioactive solution. 250 ml of the necessary eluting solution is supplied with each generator and further quantities can be ordered upon request.

The Isotope Generator contains up to 370 kBq Cs-137 as the long-lived parent nuclide which has a half life of 30.25 years. Cs-137 decays, by the emission of beta radiation, into the stable isotope Ba-137. This transition is partly effected (approx 5 %) by direct conversion into the stable energy state of Ba-137 and partly (approx 95 %) via the metastable energy state of Ba-137m.

Accompanied by the emission of gamma radiation ($E_{\gamma} = 662 \text{ keV}$) the Ba-137m decays with a half life of only 2.6-minutes into the stable isotope Ba-137.

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Isotope Generator Cs-137/Ba-137m

Handling Instructions

During the elution the Ba-137m is "milked" out of the Isotope Generator. During a one-hour educational lesson the Isotope Generator can be eluted several times. The Ba-137m activity of the eluate however, is dependent upon the radioactive equilibrium between the parent nuclide Cs-137 and the daughter nuclide Ba-137m. The maximum Ba-137m activity can only be achieved if radioactive equilibrium is reached. After each elution some time is required before the radioactive equilibrium is recovered.

After elution, the re-establishment of equilibrium can be indicated as follows:

after	1.0 min	approx 23 %
after	2.6 min	approx 50 %
after	5.2 min	approx 75 %
after	12.0 min	96 %
after	15.0 min	98 %
after	20.0 min	99.5 %

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After elution the recovery of the radioactive equilibrium within the generator can be detected with a suitable measurement instrument.

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The generator can also be used for experimental purposes as a gamma emitter. The beta radiation emitted by the Cs-137 is not detectable outside the generator capsule.

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Handling Instructions

■ Packaging

The Isotope Generator is supplied in a box containing the generator with protecting end-caps, one syringe and one plastic tube.

This box has a printed radiation warning label and can be used as the primary storage container.

The eluting solution is supplied in a separate plastic or glass bottle.

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■ Eluting procedure*

- pour a sample of the eluate into a stable flat-bottomed glass flask or into a test-tube fixed by burette pincers in a suitable stand.

- connect the plastic tube onto the syringe, draw approx 2 ml of the eluting solution into the syringe and remove the tube.

- remove the protecting cap from the generator (blue-labelled screw socket) and connect the generator with the syringe tip by tightening it lightly into the screw socket.

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Isotope Generator Cs-137/Ba-137m

Handling Instructions

- to start the eluting process, pull the protecting cap off the discharging socket (no screw thread!) and hold the generator with its discharging hole facing downwards above the receiving vessel. By carefully pressing the syringe piston the eluting solution is forced through the generator and discharged into the receiving vessel. The piston should not be pressed too strongly. The whole eluting procedure should take approximately 10 - 20 seconds.
- after the elution process has been completed, the generator should be sealed with both sealing caps. The generator should then be returned to its storage container if prolonged periods of non-usage are envisaged.
- * The colourless eluting solution has been tinted yellow for photographic purposes.

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Handling Instructions

■ Working Life

The working life of the Isotope Generator is limited by impurities which might obstruct the filters after some time, usually getting into the generator with the eluting solution. Using the eluting solution supplied by Eckert & Ziegler Nuclitec, the generator can be eluted up to 1,000 times.

■ Handling Instructions

Handling of the Eluate

Approximately 30 minutes after the elution the residual activity of the Ba-137m is lower than a thousandth of its initial activity. The total Cs-137 impurities in the eluate is < 50 Bq/ml.

The eluate needs to be disposed off according to national regulations.

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In other countries similar regulations will apply. Please consult your national Competent Authority for details.

After completion of the radioactive experiments, in compliance with the StrlSchV, the generator has to be stored in a burglar-proof cupboard which can be locked and is used for no other purpose.

If you are in doubt over any of the above handling procedures, contact your local Competent Authority for advice or call your local Eckert & Ziegler Nuclitec representative.

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Documents Supplied

The delivery of each generator includes the Handling Instructions and the Test Report. These documents must be looked after carefully.

Disposal

The Isotope Generator has to be disposed of in compliance with local legal regulations.

Isotopengenerator
Cs-137/Ba-137m

Isotope Generator
Cs-137/Ba-137m

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Handling Instructions

© Eckert & Ziegler Nuclitec
GmbH 2011 - All rights reserved.

All goods and services are sold
subject to the terms and condi-
tions of sale of Eckert & Ziegler
Nuclitec GmbH which supplies
them. A copy of these terms and
conditions is available on re-
quest.

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