


the wall thickness
of the beaker is 2mm

the filling height H depends
on the the activate volume

Scale 1:2	Tolerance ISO 2768	Surface	Material polypropylene (beaker)	H.dft	wall thickness	24.04.2018	DSt
				F	firm logo	05.01.2010	DSt
				Issue	Change	Date	Name

Nuclide various nuclides	Document Verkaufszeichnung	Responsible Dept. isotrak	Drawn DStapper	Creation Date 24.04.2018	Approval <i>[Signature]</i>
 Eckert & Ziegler		GAMMA REFERENCE SOURCE 1l Marinelli beaker Type 1N			VZ-1262-001
					Rev. H

EDV No.: \\aeant6.aea-germany.local\data2\QSAProduction\public\exchange\Zeichnungen\1251-1500\VerkaufsZeichnung\VZ-1262-001H.dft

The contents of this drawing and its enclosures are our property. The drawing and its enclosures may not be duplicated without our written approval not be made accessible to any third party. And unauthorised usage is unlawful and will be prosecuted.