


homogeneously activated
density: approx. 0,7-3g/cm³



the wall thickness
of the beaker is 2mm

the filling height H depends
on the active volume

Scale	Tolerance	Surface	Material	l.dft	wall thickness	24.04.2018	DSt
1:2	ISO 2768-		polypropylene (beaker)	G	firm logo	18.06.2009	DSt
				Issue	Change	Date	Name

Nuclide	Document	Responsible Dept.	Drawn	Creation Date	Approval	
various nuclides	Verkaufszeichnung	isotrak	DStapper	24.04.2018		
		GAMMA REFERENCE SOURCE 1l Marinelli beaker Type 1W			VZ-1520-001	
					Rev.	Issue Date
				1	24.04.18	1/1

EDV No.: \\aeant6.aea-germany.local\data2\QSAProduction\public\exchange\Zeichnungen\1501-1750\VerkaufsZeichnung\VZ-1520-001.dft

The contents of this drawing and its enclosures are our property. The drawing and its enclosures may not be duplicated without our written approval not be made accessible to any third party. And unauthorised usage is unlawful and will be prosecuted.